

BPC-150

BI-POLAR ELECTRO CAUTERY



BPC-150 BI-POLAR CAUTERY

- High frequency solid state circuitry.
- Compact in size, elegant look and easy to use.
- Reusable forceps, erasers and cables.
- Audible indication.
- Dial for 0-9 Variable output power control by nab.
- Control by foot switch with buzzer sound.
- Short circuit protection.
- Autoclavable Silicon cable with Forcep/Eraser
- Imported Plastic box.

OPTIONAL ACCESSORIES



LONG CABLE



ERASERS



CURVED COATED
BI-POLAR FORCEP



STRAIGHT COATED
BI-POLAR FORCEP

TECHNICAL SPECIFICATION

Supply Voltage:	230V AC or 110V AC
Frequency:	230V AC or 110V AC
HF Output Power:	10Watts
Power Adjustment Control:	Power Adjustment Control
Display Range:	0-20
Dimension:	140(W) x 50(H) x 170(L)
Casing:	ABS plastic / M.S powder coated
Weight:	1.7 Kg
Size:	225mmX200mmX72mm
Activation of HF Power:	Foot Switch(detachable)



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008



www.mediequipbazar.com