



# A breath of fresh air...

## ...the new slimline Pulsair Desktop Tonometer

A new dimension in glaucoma  
screening elegance

- Small and space saving footprint combined with the elegant and slim optical mainframe allows the clinician to screen all patients with comfort and respect
- Fast, accurate and easy to use
- Slim and open design with clear user controls
- Soft and comfortable puff



**Keeler**  
*Further details overleaf ►*

# The new slimline Pulsair Desktop Tonometer

For over 90 years Keeler has continued to inn-  
design and manufacture world class ophthaln  
instruments in our Windsor factory.

In 1986 the Pulsair Tonometer revolutionised  
contact tonometry with its unique hand held  
and now the new Pulsair Desktop brings the l  
name to those requiring a chinrest desktop d

The Pulsair Desktop has a small and space sav  
footprint that, combined with the elegant an  
optical mainframe, allows it to blend seamles  
the clinical environment. The openness of the  
increases the confidence of both patient and

It is also uncomplicated and therefore quick t  
the novice and professional alike. Taking cont  
tonometry has never been easier or faster. Cle  
controls and a colour video alignment screen  
combine to create a thoroughbred tonometer

Contact us today for further details.



## Specifications

<b>Console dimensions</b>	475 x 410 x 245mm (H x D x W)
<b>Console weight</b>	16kg
<b>Calibrated range</b>	5mmHg to 50mmHg
<b>Working distance</b>	20mm from surface of patient's cornea to front surface of first lens. This equates to a nominal distance of 15mm from the front of the puff tube shroud to the front surface of the patient's cornea
<b>Displayed scale</b>	Single line 16 character alphanumeric display
<b>Illumination system</b>	LED infra red
<b>Complies with</b>	Electrical Safety (Medical) BS EN 60601-1:2006 Electromagnetic compatibility EN 60601-1-2:2007 Ophthalmic instruments - fundamental requirements and test methods ISO 15004-1:2006 Optical radiation hazard ISO 15004-2:2007
<b>Power Supply Unit</b>	Switch mode, (110 – 240V)+/- 10% multi plug type compliant to EN 60601-1 EN 61000-6-2, EN 61000-6-3
<b>Power supply output</b>	30 VA (12V DC 2.5A)
<b>Frequency</b>	50/60 Hz
<b>Warranty</b>	2 years

Authorized Sales Partner:

#### HEAD OFFICE

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Plot No. 50/1B-50/2B, Rajharsh  
Colony, Near CI Square Shopping  
Mall, Kolar Road, Bhopal, M.P  
India-462042

#### DELHI SALES OFFICE

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Building Name - Classic's, T-1895/1,  
Upper Ridge Road, Karol Bagh, New  
Delhi India- 110005

#### RAIPUR SALES OFFICE

**MEX ENTERPRISES**  
ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
House No.-10/806, ShivanandNagar,  
Jhanda Chowk, Raipur Chhattisgarh  
India-492008