

# PRIMA OPH



## GIVING SURGEONS THE GIFT OF LIGHT



With increasing demands on ophthalmic surgeons to provide better and more affordable care, they must be empowered with instruments that perform with ultimate precision. The Prima OPH is the first in a new class of microscopes designed for doctors with a focus on better and more affordable eye care.


With LABOMED's proprietary *MaxLite™* coated optics, superb images with true-to-life rendition and glare-free view are guaranteed. A sturdy, balanced stand and a host of ergonomic design features deliver ultimate user convenience.

## ULTRA BRIGHT LED + LONG-LASTING BULB = DELIGHT

**1** Bright and even illumination coupled with apochromatic optics and an ergonomic stand makes the Prima OPH an ideal ophthalmic solution.



**2** Red reflex is made possible in the Prima OPH thanks to an ingenious optical design enabling visualization of an illuminated retina for easy detection of ruptured vessels, arteries, and damage to the optic disc.



**3** LABOMED developed the *ProLine* range of adapters for DSLRs and video cameras, enabling convenient digital documentation and live viewing. The *iVu* series is available to doctors who want a camera that does it all: live HD streaming, image and video capture, and imaging software. For assistant viewing, the LABOMED assistoscope is available with an ergonomic position and light-weight construction for a well-balanced system.



# PRIMA OPH

## SPECIFICATIONS

### Binocular tubes Eyepieces

### Achromatic magnifier Objective

### Fine focus

### Light source

### Built-in filters

### Vertical movement of arm

### Microscope carriers

### Optional accessories

### PRIMA OPH X-Y & Z AXIS

0-210° tiltable ergonomic head  
WF 10X/16 mm with eye guards,  
diopter adjustment  $\pm 5$  mm and  
diopter lock

Optional: WF 12.5X/16 mm

5 steps: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X  
f-200 mm

Optional: f-175 mm

Range 40 mm, motorized

50W LED

Green and Blue

550 mm

Multi-axial (XYZ) foot control

45° inclined head

Double beam splitter

Assistant binocular attachment

*ProLine* range of camera adapters

*iVu* S5 5MP integrated digital  
camera module

### PRIMA OPH Z AXIS

0-210° tiltable ergonomic head  
WF 10X/16 mm with eye guards,  
diopter adjustment  $\pm 5$  mm and

Optional: WF 12.5X/16 mm

5 steps: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X  
f-200 mm

Optional: f-175 mm

Range 40 mm, motorized

50W LED

Green and Blue

550 mm

Single axis (Z) foot control


45° inclined head

Double beam splitter

Assistant binocular attachment

*ProLine* range of camera adapters

*iVu* S5 5MP integrated digital  
camera module



### DESCRIPTION

Prima OPH

Prima OPH

### CONFIGURATION

Basic, motorized fine focus

Motorized fine focus and  
X-Y coupling

### CATALOG NO.

6136000-000

6135000-000



**LABOMED**

920 Auburn Aourt  
Fremont, CA 94538



Authorized Sales Partner:

**MEX**  
**ENTERPRISES**  
*We Share in Healthcare...*

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

### HEAD OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Plot No. 50/1B-50/2B, Rajharsh  
Colony, Near CI Square Shopping  
Mall, Kolar Road, Bhopal, M.P  
India-462042

### DELHI SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
Building Name - Classic's, T-1895/1,  
Upper Ridge Road, Karol Bagh, New  
Delhi India- 110005

### RAIPUR SALES OFFICE

#### MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company  
House No.-10/804, ShivanandNagar,  
Jhanda Chowk, Raipur Chhattisgarh  
India-492008

Ph.: +91 755 4256391, +91 9302778762, +91 9111429222, Email: info@mediequipbazar.com, export@mediequipbazar.com

Website: www.mediequipbazar.com