

MEL-900

EXAMINATION LIGHT



(Wall Mounted)



(Wall Mounted with Extended Rod)



Optional Battery Backup



FEATURES



- **Shadowless Light.**
- **Touch Controller.**
- **With special European design ventilation of heat.**
- **Body is made of high quality resistant polymer.**

TECHNICAL SPECIFICATION

Intnsity	55,000 Lux+/-10%
Size of Light Field	125 cm
Colour Temperature	4500-5000k
Colour Rendering Index	93 RA
LED Life	>50,000 Hrs.
Number of LED	21
Diameter of Light	150 mm
Brightness Control	Touch
Power Supply	Power 220 Volts/50 Hz AX
Battery Backup	Optional
Colour of LED	Double(Yellow & White)
Illumination adjustment	0% to 100%



+91 9302 778 762
 +91 9111 429 222



export@mediequipbazar.com
 info@mediequipbazar.com



+91 9302 778 762
 +91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
 Industrial Park Phase I,
 Badiyakhedi, Sehore, M.P.,
 India - 466 001

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
 Near CI Square Shopping Mall, Kolar Road,
 Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
 Upper Ridge Road, Karol Bagh,
 New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
 Jhanda Chowk, Raipur CG-492008



MANUFACTURING UNIT

Plot No.-40, KI Ansal
 Green, Kolar Road
 Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
 Shade, Near D-Mart,
 Bhopal, 462042



www.mediequipbazar.com

