

MEDIX
MEDICAL

MEX
ENTERPRISES
We Share in Healthcare...

HCL-5000

PORTABLE SLIT LAMP




www.mediequipbazar.com

HCL-5000 HAND HELD SLIT LAMP

Medix HCL-5000 The Portable Slit Lamp utilizes a high luminance white LED light source, coupled with a lightweight, compact design, which could give you the comfortable of using a slit lamp in the palm of your hand. It is the best choice for non-ambulatory patient visits, off-site clinics, examinations of young children and patients who don't fit into the traditional examination lane.

- High Illumination White LED Light Source
- Rechargeable Green Lithium Battery
- Fully Charged Can Be Used for 5 or 6 Hours
- Light-weight & Compact Design, Single Hand Operate
- Link-able Eyepiece Equip With Four Spots
- Can Be Digitized By Coopix Camera or Phone Eyepiece Adapters
- Packing Include: 1*16X/10X Eyepiece; 1*Forehead Rest; 1*Focusing Stick; 1*Lithium Battery; 1*Battery Charger; 1*Packing Box





TECHNICAL SPECIFICATION

Type:	Binocular Crossed Stereo Microscope
Eyepiece:	10X, 16X
Total Magnification:	10X, 16X(Optional)
Diopter Compensation:	+/-5D
Slit Width:	0~10mm(Adjustable)
Spot Diameter:	0.5mm, 3mm, 5mm, 10mm
Filter:	Head Absorption, Grey, Red-Free, Cobalt Blue
Working Time:	5~6 Hours(Constantly)
Working Distance:	80mm
Net Weight:	750g
Gross Weight:	2600g



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008

