

MEDIX MED-200

ENT SURGICAL MICROSCOPE



The MED-200 microscope manufactured by Mex Enterprises fulfils the most demanding requirements for ENT examination and diagnosis.

Its excellent optical and mechanical features, together with its great mobility and versatility, make the MED-200 microscope the ideal equipment for any treatment.

TECHNICAL SPECIFICATION

Technical Specification

Working distance: F=250mm objective lens(400mm-optional)

Eye pices: 12.5x wide field

Optical head: Straight viewing type 90 degree

PD Adjustment: 55mm to 75mm

Magnification: 0.4x,0.6x,1.0,1.6x,2.5x (Five step)

Magnification: 0.6x,1.0,1.6x, (Three step)

Illumination Section

Light Source: Halogen 12V 100W (80000 lux)

Filter: Build in heat absorbing filter and UV,
green and blue filters with switchable facility

Illumination Control: Continuous Adjustment

Stand Section

Type: Floor Stand

Rotation of Arm: 320° with lock

Vertical Height Adjustment



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No.-79,
Commercial Shade,
Near D-Mart,
Bhopal-462042 MP



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No.-59, KI Estate,
Ahmedabad-382430
Gujrat



www.mediequipbazar.com

REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near Cl Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008

