

MED-200

OPHTHALMIC SURGICAL MICROSCOPE



IMPORTED OPTIC HEAD

 **MADE
IN
GERMANY**

■ MED-200

The MED-200 microscope manufactured by Medix Medical fulfils the most demanding requirements for Ophthalmic Surgery.

Its excellent optical and mechanical features, together with its great mobility and versatility, make the MED-200 microscope the ideal equipment for any treatment.

■ DESIGN

A bright Halogen light source delivers a natural daylight color temperature. With Standby spare illumination slot provides quick swap for light source as needed. Compact enclosure for electronics for easy movement.



■ OPTICS

An German made imported apochromatic optical system produces superior resolution and artifact-free images for excellent visibility. Optics coatings promise reflection-free, durable optics.

■ SUPPORTS

Rolling Floor, Stand Fixed Floor,
Mount Wall, Mount Ceiling Mount

MEDIGLOW

We offer this microscope with a detachable Mediglow, a brilliant red glow system which provides a omnipresent Red glow (Retro illumination) during Surgery which makes Capsulorhexis more visible and also increase depth of field.

It comes with separate light source so it could be used independently from microscope.



TECHNICAL SPECIFICATION

Eyepieces

45° inclined head
WF 10X/16 mm with eye guards,
diopter adjustment ± 5 mm and
diopter lock

Apochromatic magnichanger Objective

Optional: WF 12.5X/16 mm
5 steps: 0.4X, 0.6X, 1.0X, 1.6X, 2.5X
f-200 mm

Fine focus

Optional: f-175 mm

Light source

Range 40 mm, motorized

Built-in filters

12V 100W Halogen Lamp

Vertical movement of arm

Green and Blue

Microscope carriers

550 mm


Optional accessories

Multi-axial (XYZ) foot control

Double beam splitter

Assistant binocular attachment

range of camera adapters



MANUFACTURED BY

MEX
ENTERPRISES
We Share in Healthcare...

HEAD OFFICE MEX ENTERPRISES

Plot No. 50/1B-50/2B, Rajharsh
Colony, Near CI Square Shopping
Mall, Kolar Road, Bhopal, M.P
India-462042

DELHI SALES OFFICE MEX ENTERPRISES

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh, New
Delhi India- 110005

RAIPUR SALES OFFICE MEX ENTERPRISES

House No.-10/806, ShivanandNagar,
Jhanda Chowk, Raipur Chhattisgarh
India-492008

Ph.: +91 755 4256391, +91 9302778762, +91 9111429222,

Email: info@mediequipbazar.com, export@mediequipbazar.com Website: www.mediequipbazar.com