

PRISM BAR

MODEL:VP-15 & HP-16



PRISM BAR

Optel high quality acrylic prism available in prism bar the doctor looks for any movement in the eye that is uncovered.

The prism and alternate cover test alternate the cover over both eyes while a prism is placed in front of one eye.

This measures the offset(difference) between the two eyes to find the right prism to fix your double vision.

PRISM BAR

Product Parameters

Model: VP-15 & HP-16

- Two pos in one aluminium case
- Vertical Diopter: 1 2 3 4 5 6 7 8 9 10 12 14 16 18 20 25 30
- Horizontal Diopter: 1 2 4 6 8 10 12 14 16 18 20 25 30 35 40 45
- Material: Acrylic
- Aluminium case



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008