


# R.O.S.

## RETINA OBSERVATION SYSTEM



## FEATURES

- Integrated inverter which minimize the height of the microscope and flips away when needed
- Ultra view non inverter optical system
- Easy to adapt with any microscope
- Excellent field of view with better stereoscopic vision
- High optical resolution with focussing mechanism to easy the focus of the retina
- Optimum rotation without retinal shift
- Uninterrupted access to surgeon at higher diopter



## TECHNICAL SPECIFICATION

<b>Microscope Working Distance:</b>	175mm or 200mm
<b>Front Mid Field Lens:</b>	for approx 90 degree viewing angle
<b>Front Wide Angle Lens:</b>	for approx 124 degree viewing angle



+91 9302 778 762  
 +91 9111 429 222



export@mediequipbazar.com  
 info@mediequipbazar.com



+91 9302 778 762  
 +91 9111 429 222



**MANUFACTURING UNIT**

Plot No. 23-B, Badiyakhedi  
 Industrial Park Phase I,  
 Badiyakhedi, Sehore, M.P.,  
 India - 466 001

**REGISTERED OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,  
 Near Cl Square Shopping Mall, Kolar Road,  
 Bhopal, M.P.-462042

**SALES OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,  
 Upper Ridge Road, Karol Bagh,  
 New Delhi - 110005

**SALES OFFICE**

**MEX ENTERPRISES**

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,  
 Jhanda Chowk, Raipur CG-492008



**MANUFACTURING UNIT**

Plot No.-40, KI Ansal  
 Green, Kolar Road  
 Bhopal-462042



**MANUFACTURING UNIT**

Plot No. - 79, Commercial  
 Shade, Near D-Mart,  
 Bhopal, 462042



www.mediequipbazar.com

