



VERTICAL AUTOCLAVES & DRESSING DRUM

300mm X 400mm + 150mm Lid. (12"X22")-2KW Load--39 Ltrs.With Lid



VERTICAL AUTOCLAVE PRODUCT DETAILS

Introducing the Medix Electric Autoclave Pressure Aluminium Body, a highly reliable and efficient solution for sterilization. As one of the largest and most trusted manufacturers, suppliers and exporters of Autoclaves in the world. Medix Medical takes pride in delivering top-quality Autoclaves to hospitals, healthcare facilities, clinics, food processing units, tissue culture labs and laboratories worldwide. With our commitment to excellence, we supply Autoclaves to numerous countries across the globe. Let's explore the exceptional features and specifications of our Autoclave:

- **Design:** Our Autoclave features a pressure cooker design that allows for ease of use. This ensures convenience and efficiency during the sterilization process, making it suitable for various applications.
- **Hi-Quality Aluminium Construction:** All our Autoclaves are constructed from heavy-gauge aluminium, guaranteeing durability and reliability. The seamless aluminium construction provides a bacteria-free environment, ensuring optimal sterilization.
- **Sterilization Cycle:** Designed for a sterilizing cycle of 25-30 minutes at a temperature range of 120°C to 135°C, our Autoclave delivers effective sterilization within a short duration. This helps streamline workflow and improve product in medical and laboratory settings.
- **Advanced Safety Features:** Our Autoclave is equipped with a two-stage over-pressure protection system, incorporating calibrated continuous bleeding pressure stopcock and a pressure release weight valve. These features ensure reliable secure operations, providing peace of mind.
- **Precise Monitoring:** The Autoclave is equipped with a color-coded pressure gauge, allowing easy monitoring of the internal chamber conditions during the sterilization cycle. The sterilization zone is highlighted in green, indicating the pressure range of 15 psi to 20 psi (or 1.034 bar to 1.378 bar), equivalent to a temperature range of 121°C to 127°C. The caution zone is marked in red for quick visual reference.
- **Convenient and Safe Handling:** Our Autoclave features heat-resistant handles and a steam release/control valve, ensuring safe and convenient handling during the sterilization process. These design elements prioritize user safety and comfort.

STAINLESS STEEL DRESSING DRUM PRODUCT DETAILS

DRESSING DRUM

Medix is a reputable manufacturer, supplier, and exporter of Dressing Drum designed for the effective sterilization and storage of medical instruments and equipment. These Dressing Drums made from high-quality stainless steel and crafted with heavy gauge stainless steel sheets, these drums ensure durability and longevity.

- The Dressing Drums are available in 11"X9" sizes, including. They are specially designed to fit inside autoclaves, allowing for convenient sterilization processes. With their robust construction and reliable performance, these Dressing Drums provide a safe and hygienic storage solution for medical instruments, ensuring optimal sterilization and ease of use.
- Choose Medix Medical for high-quality Dressing Drums that international standards, offering exceptional performance and longevity for medical facilities, hospitals and healthcare professionals.



+91 9302 778 762
+91 9111 429 222



export@mediequipbazar.com
info@mediequipbazar.com



+91 9302 778 762
+91 9111 429 222



MANUFACTURING UNIT

Plot No. 23-B, Badiyakhedi
Industrial Park Phase I,
Badiyakhedi, Sehore, M.P.,
India - 466 001



MANUFACTURING UNIT

Plot No.-40, KI Ansal
Green, Kolar Road
Bhopal-462042



MANUFACTURING UNIT

Plot No. - 79, Commercial
Shade, Near D-Mart,
Bhopal, 462042



www.mediequipbazar.com



REGISTERED OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Plot No. 50/1B-50/2B, Rajharsh Colony,
Near CI Square Shopping Mall, Kolar Road,
Bhopal, M.P.-462042

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

Building Name - Classic's, T-1895/1,
Upper Ridge Road, Karol Bagh,
New Delhi - 110005

SALES OFFICE

MEX ENTERPRISES

ISO 13485 :2016 & ISO 9001 :2015 Certified Company

House No.-10/806, Shivanand Nagar,
Jhanda Chowk, Raipur CG-492008